

# Cell & Tissue DNA Extraction Kit

Catalog Number: EXNA075

Size: 50 Tests

For Research Use Only. Not Intended for Diagnostic or Therapeutic Use.



info@exreprotein.com

## Product Description

The Cell & Tissue DNA Extraction Kit uses a centrifugal adsorption column that can specifically bind DNA and a unique buffer system for rapid extraction of genomic DNA from animal tissues and cells. The optimized Binding Buffer and Proteinase K rapidly lyses cells and inactivates intracellular nucleases, and genomic DNA is selectively adsorbed on the silicon matrix membrane in the centrifuge column. A series of rapid rinse-centrifugation steps remove cell metabolites, proteins, and other impurities. The purified genomic DNA is compatible with Southern-blot and enzyme digestion reactions.

The kit is efficient and simple-to-use. Single samples generate highly pure genomic DNA in under 30 minutes.

## Limitations

- For Research Use Only. Not Intended for Diagnostic or Therapeutic Use.

## Precautions

- Always wear appropriate protective clothing and follow safe laboratory procedures.
- Minimize the length of UV exposure when excising the band from the gel to avoid damage to DNA.

## Materials Provided & Storage

Shipped at ambient temperature. Proteinase K shipped on blue ice.

Store the kit at 2-30°C for up to 24 months. Proteinase K in glycerol buffer, store at 2-8°C for 12 months or -15 to -25°C for 24 months.

| Component                  | Size     | Quantity |
|----------------------------|----------|----------|
| Balanced Buffer            | 5 mL     | 1 Bottle |
| Lysate Buffer              | 11 mL    | 1 Bottle |
| Binding Buffer             | 11 mL    | 1 Bottle |
| Inhibitor Removal Solution | 25 mL    | 1 Bottle |
| Wash Buffer                | 15 mL    | 1 Bottle |
| Elution Buffer             | 15 mL    | 1 Bottle |
| Proteinase K               | 1 mL     | 1 Bottle |
| Adsorption Column          | 50 Tubes | 1 Pack   |
| Collection Tube            | 50 Tubes | 1 Pack   |

## Additional Materials Required, Not Supplied

|                   |                         |
|-------------------|-------------------------|
| 1X PBS            | RNase A                 |
| Isopropanol       | 1.5 mL centrifuge tubes |
| Ethanol (96-100%) |                         |

## Preparation of Reagents

**Wash Buffer:** Add 60 mL ethanol (96-100%) to the Wash Buffer before use.

## Protocol

**Column equilibration:** The adsorption column should be pretreated with Balanced Buffer before the experiment to improve the nucleic acid binding ability. Place the Adsorption Column into the Collection Tube, add 100  $\mu$ L of Balanced Buffer to each tube pairing, centrifuged at 13,000 $\times$ g for 60 seconds, discard the waste liquid from the Collection Tube, and place the Adsorption Column back into the Collection Tube for use later.

### Cultured Cell Preparation:

1. Add  $10^6$ - $10^7$  suspended cells to 1.5 mL centrifuge tube (do not exceed  $10^7$  cells). For adherent cells, treat with Trypsin.
2. Centrifuge at 13,000  $\times$  g for 60 seconds to remove the supernatant.
3. Resuspend cells in 200  $\mu$ L of 1X PBS, centrifuge at 13,000  $\times$  g for 60 seconds, discard the supernatant, and resuspend the cells in 180  $\mu$ L of 1X PBS.
4. Add 20  $\mu$ L of Proteinase K solution, mix thoroughly, then add 200  $\mu$ L Binding Buffer, vortex thoroughly, and place at 70°C for 10 minutes. **Optional:** if RNA removal is required, add 20  $\mu$ L RNase A (25 mg/mL), shake to mix, and leave at room temperature for 5-10 minutes after the completion of step 4.
5. Proceed to Extraction Protocol.

### Animal Tissue Preparation:

1. Add 180  $\mu$ L of Lysate Buffer to fresh or thawed tissue ( $\leq$  25 mg), homogenize on ice.
2. Add 20  $\mu$ L of Proteinase K solution and vortex thoroughly.
3. Place the lysates in a water bath at 56°C for 1-3 hours or until tissue digestion is complete, shaking gently several times to assist lysis. **Optional:** if RNA removal is required, add 20  $\mu$ L RNase A (25 mg/mL), shake to mix, and leave at room temperature for 5-10 minutes after the completion of step 3.
4. Add 200  $\mu$ L Binding Buffer, vortex thoroughly, and place at 70°C for 10 minutes.
5. Proceed to Extraction Protocol.

### Extraction Protocol:

1. Allow sample to cool to room temperature. Add 100  $\mu$ L isopropanol, immediately vortex to mix thoroughly.
2. Place the Adsorption Column into the Collection Tube. Add the mixture (including the flocculent precipitate) to the Adsorption Column, centrifuge at 13,000  $\times$  g for 60 seconds, and discard the waste liquid.
3. Add 500  $\mu$ L of Inhibitor Removal Solution, centrifuge at 12,000  $\times$  g for 60 seconds, and discard the waste liquid.
4. Add 600  $\mu$ L of Wash Buffer (per Preparation of Reagents), centrifuge at 12,000  $\times$  g for 60 seconds, discard the waste liquid, and repeat this step once.
5. Place the Adsorption Column back into the Collection Tube and centrifuge at 13,000  $\times$  g for 2 minutes to thoroughly remove the Wash Buffer and minimize residual ethanol affecting the downstream reaction.
6. Place the Adsorption Column in a new 1.5 mL centrifuge tube, add 100  $\mu$ L Elution Buffer to the center of the membrane and place it at room temperature for 3-5 minutes.
7. Elute genomic DNA by centrifugation at 12,000  $\times$  g for 60 seconds and store at -20°C. **Note:** the eluted solution can be pipetted back into the Adsorption Column, allowed to stand for 3-5 minutes, and centrifuged at 12000  $\times$  g for 60 seconds to increase the DNA concentration. Preheating the Elution Buffer to 65-70°C in advance may improve results.